

# Clinical Accuracy of TRITEMP™ (non-contact infrared forehead thermometer)

EN ISO 80601-2-56:2017 and ASTM E1965-98:2016 are both FDA-recognised voluntary consensus standards used by device manufacturers to evaluate the performance of non-contact infrared thermometers by (i) testing the accuracy of the device against a standard blackbody source and (ii) performing a clinical study to evaluate the accuracy and the effectiveness of the device in clinical settings (*Sullivan et al., 2021*).

The CLINICAL ACCURACY of TRITEMP™ with REFERENCE CLINICAL THERMOMETERS widely used in hospitals was performed in accordance with the requirements of both EN ISO 80601-2-56 and ASTM E1965-98:2016.

## TEST RESULTS

The results of CLINICAL ACCURACY VALIDATION are summarised in table below including the CLINICAL BIAS ( $\Delta_{cb}$ ), with its LIMITS OF AGREEMENT (LA) and the CLINICAL REPEATABILITY ( $\sigma_r$ ).

	ISO 80601-2-56	ASTM E1965-98	
	All (n=771)	Paediatrics (n=253)	Adults (n=518)
CLINICAL BIAS ( $\Delta_{cb}$ )	0.08°C	0.03°C	0.10°C
LIMITS OF AGREEMENT (LA)	0.71°C	0.55°C	0.78°C
CLINICAL REPEATABILITY ( $\sigma_r$ )	0.07°C	0.09°C	0.06°C

Note: The reference body site was axilla and oral for paediatrics and adults, respectively.

**CLINICAL BIAS ( $\Delta_{cb}$ )** – Mean difference between output temperatures of a clinical thermometer and a reference clinical thermometer for the intended reference body site with specified limits of agreement when measured from selected group of subjects.

**LIMITS OF AGREEMENT (LA)** – The magnitude of a potential disagreement between outputs of two clinical thermometers equal to double the standard deviation of output temperature differences when used on the same human subject.

**CLINICAL REPEATABILITY ( $\sigma_r$ )** – Pooled standard deviation (over a selected group of subjects) of changes in multiple output temperatures taken from the same subject at the same measuring site with the same clinical thermometer by the same operator within a relatively short time.

## CONCLUSION

From the above results, it can be concluded that this product, a non-contact infrared forehead thermometer TRITEMP™, can be used for measuring body temperatures accurately within a clinical setting.

Clinical accuracy was validated in accordance with ISO 80601-2-56 on 771 subjects, adults and paediatrics, against a reference clinical thermometer for oral (adults) and axilla (paediatric) reference sites. Results found a clinical bias of 0.08°C with limits of agreement of 0.71°C and a clinical repeatability of 0.07°C.

Clinical accuracy was validated in accordance with ASTM E1965-98 on 771 subjects, adults and paediatrics, against a reference clinical thermometer for oral (adults) and axilla (paediatric) reference sites. Results found a clinical bias of 0.10°C with limits of agreement of 0.78°C and a clinical repeatability of 0.06°C for 518 oral (adult) measurements and a clinical bias of 0.03°C with limits of agreement of 0.55°C and a clinical repeatability of 0.09°C for 253 axilla (paediatric) measurements.